
Compatibility Chart



■ compatible x not compatible

ahn
BVS
www.bvskorea.com
+82-02-546-6212
+82-02-546-6211

Compatibility Chart

AHN — Pipet4u Pro / Pipet4u Revolution — single and multi channel pipettes

Description	Cat. No.	single 0.2-2 µL	single 0.5-10 µL	single 2-20 µL	single 5-50 µL	single 10-100 µL	single 20-200 µL	single 100-1000 µL	single 1000-5000 µL	single 1000-10000 µL	multi 0.5-1.0 µL	multi 5-50 µL	multi 10-100 µL	multi 20-300 µL
AHN myTip® Filter Tips														
AHN myTip® FT 10 µL, clear	2-002-81-0	■	■								■			
AHN myTip® FT 10 µL, clear, sterile	2-003-C5-0	■	■								■			
AHN myTip® FT 10/20 µL XL, clear	2-062-81-0	■	■	■							■			
AHN myTip® FT 10/20 µL XL, clear, sterile	2-063-C5-0	■	■	■							■			
AHN myTip® FT 20 µL, clear	2-117-80-0				■	■	■					■	■	
AHN myTip® FT 20 µL, clear, sterile	2-118-C5-0				■	■	■					■	■	
AHN myTip® FT 100 µL, clear	2-119-80-0					■	■					■	■	
AHN myTip® FT 100 µL, clear, sterile	2-120-C5-0					■	■					■	■	
AHN myTip® FT 200 µL, clear	2-127-80-0					■	■					■	■	
AHN myTip® FT 200 µL, clear, sterile	2-128-C5-0					■	■					■	■	
AHN myTip® FT 300 µL, clear	2-134-80-0					■	■					■	■	■
AHN myTip® FT 300 µL, clear, sterile	2-135-C5-0					■	■					■	■	■
AHN myTip® FT 1000 µL, clear	2-202-80-0							■						
AHN myTip® FT 1000 µL, clear, sterile	2-203-C4-0							■						
AHN myTip® FT 1250 µL, clear	2-204-80-0							■						
AHN myTip® FT 1250 µL, clear, sterile	2-205-C4-0							■						
AHN myTip® LT Low Retention Pipette Tips														
AHN myTip® LT 10 µL, clear	4-002-81-0	■	■								■			
AHN myTip® LT 10 µL, clear, sterile	4-003-C5-0	■	■								■			
AHN myTip® LT 10/20 µL XL, clear	4-062-81-0	■	■	■							■			
AHN myTip® LT 10/20 µL XL, clear, sterile	4-063-C5-0	■	■	■							■			
AHN myTip® LT 200 µL, clear	4-121-86-0				■	■	■					■	■	
AHN myTip® LT 200 µL, clear, sterile	4-122-C5-0				■	■	■					■	■	
AHN myTip® LT 300 µL, clear	4-127-80-0				■	■	■					■	■	■
AHN myTip® LT 300 µL, clear, sterile	4-128-C5-0				■	■	■					■	■	■
AHN myTip® LT 1000 µL, clear	4-202-80-0							■						
AHN myTip® LT 1000 µL, clear, sterile	4-203-C4-0							■						
AHN myTip® LT 1250 µL, clear	4-204-80-0							■						
AHN myTip® LT 1250 µL, clear, sterile	4-205-C4-0							■						
AHN myTip® LFT Low Retention Filter Tips														
AHN myTip® LFT 10 µL, clear, sterile	5-003-C5-0	■	■								■			
AHN myTip® LFT 10/20 µL XL, clear, sterile	5-063-C5-0	■	■	■							■			
AHN myTip® LFT 20 µL, clear, sterile	5-118-C5-0				■	■	■					■	■	
AHN myTip® LFT 100 µL, clear, sterile	5-120-C5-0					■	■					■	■	
AHN myTip® LFT 200 µL, clear, sterile	5-128-C5-0					■	■					■	■	
AHN myTip® LFT 300 µL, clear, sterile	5-135-C5-0					■	■					■	■	■
AHN myTip® LFT 1000 µL, clear, sterile	5-203-C4-0							■						
AHN myTip® LFT 1250 µL, clear, sterile	5-205-C4-0							■						

Gilson — Pipetman / Pipetman G — single and multi channel pipettes

Cat. No.	single 0.2-2 µL	single 0.5-10 µL	single 2-20 µL	single 10-100 µL	single 20-200 µL	single 100-1000 µL	single 1000-5000 µL	single 1000-10000 µL	multi 10 µL	multi 20 µL	multi 200 µL	multi 300 µL
2-002-81-0	■	■	x						x	x		
2-003-C5-0	■	■	x						x	x		
2-062-81-0	■	■	x						x	x		
2-063-C5-0	■	■	x						x	x		
2-117-80-0			x						x	x		
2-118-C5-0			x						x	x		
2-119-80-0			x	■					x	x		
2-120-C5-0			x	■					x	x		
2-127-80-0			x	x	■				x	x	x	
2-128-C5-0			x	x	■				x	x	x	
2-134-80-0					■						■	■
2-135-C5-0					■						■	■
2-202-80-0						■						
2-203-C4-0						■						
2-204-80-0						■						
2-205-C4-0						■						
4-002-81-0	■	■	x						x	x		
4-003-C5-0	■	■	x						x	x		
4-062-81-0	■	■	x						x	x		
4-063-C5-0	■	■	x						x	x		
4-121-86-0			■	■	■				■	■	x	
4-122-C5-0			■	■	■				■	■	x	
4-127-80-0			■	■	■					■	■	■
4-128-C5-0			■	■	■					■	■	■
4-202-80-0						■						
4-203-C4-0						■						
4-204-80-0						■						
4-205-C4-0						■						
5-003-C5-0	■	■	x						x	x		
5-063-C5-0	■	■	x						x	x		
5-118-C5-0			x						x	x		
5-120-C5-0			x	■					x	x		
5-128-C5-0			x	x	■				x	x	x	
5-135-C5-0					■						■	■
5-203-C4-0						■						
5-205-C4-0						■						

Eppendorf — Research Plus — single and multi channel pipettes

Cat. No.	single 0.2-2 µL	single 0.5-10 µL	single 2-20 µL (grey)	single 2-20 µL (yellow)	single 10-100 µL	single 20-200 µL	single 100-1000 µL	single 500-5000 µL	single 1000-10000 µL	multi 10 µL	multi 100 µL	multi 300 µL
2-002-81-0	■	■								■		
2-003-C5-0	■	■								■		
2-062-81-0	■	■	■	x						■		
2-063-C5-0	■	■	■	x						■		
2-117-80-0			x	■								
2-118-C5-0			x	■								
2-119-80-0			x	■	x						■	
2-120-C5-0			x	■	x						■	
2-127-80-0			x	x	x	■					■	
2-128-C5-0			x	x	x	■					■	
2-134-80-0						■						■
2-135-C5-0						■						■
2-202-80-0							■					
2-203-C4-0							■					
2-204-80-0							■					
2-205-C4-0							■					
4-002-81-0	■	■								■		
4-003-C5-0	■	■								■		
4-062-81-0	■	■	■	x						■		
4-063-C5-0	■	■	■	x						■		
4-121-86-0			x	■	x	■					■	
4-122-C5-0			x	■	x	■					■	
4-127-80-0				■	■	■					■	■
4-128-C5-0				■	■	■					■	■
4-202-80-0							■					
4-203-C4-0							■					
4-204-80-0							■					
4-205-C4-0							■					
5-003-C5-0	■	■								■		
5-063-C5-0	■	■	■	x						■		
5-118-C5-0			x	■								
5-120-C5-0			x	■	x						■	
5-128-C5-0			x	x	x	■					■	
5-135-C5-0						■						■
5-203-C4-0							■					
5-205-C4-0							■					

Brand — Transferpette, Transferpette S — single and multi channel pipettes

Cat. No.	single 0.5-10 µL S	single 2-20 µL S	single 10-100 µL	single 20-200 µL	single 100-1000 µL	multi 10 µL S	multi 100 µL S	multi 200 µL	multi 300 µL S
2-002-81-0	■					■			
2-003-C5-0	■					■			
2-062-81-0	■					■			
2-063-C5-0	■					■			
2-117-80-0		■							
2-118-C5-0		■							
2-119-80-0			x				■		
2-120-C5-0			x				■		
2-127-80-0			x	■				■	
2-128-C5-0			x	■				■	
2-134-80-0				■				■	■
2-135-C5-0				■				■	■
2-202-80-0					■				
2-203-C4-0					■				
2-204-80-0					■				
2-205-C4-0					■				
4-002-81-0	■					■			
4-003-C5-0	■					■			
4-062-81-0	■					■			
4-063-C5-0	■					■			
4-121-86-0		■	x	■			■	■	
4-122-C5-0		■	x	■			■	■	
4-127-80-0		■	■	■			■	■	■
4-128-C5-0		■	■	■			■	■	■
4-202-80-0					■				
4-203-C4-0					■				
4-204-80-0					■				
4-205-C4-0					■				
5-003-C5-0	■					■			
5-063-C5-0	■					■			
5-118-C5-0		■							
5-120-C5-0			x				■		
5-128-C5-0			x	■				■	
5-135-C5-0				■				■	■
5-203-C4-0					■				
5-205-C4-0					■				

Thermo — Finnpipette F1 / Novus — manual single and multi channel and electronic single and multi channel pipettes

Cat. No.	manual single 0.5-10 µL pink	manual single 0.5-10 µL yellow	manual single 2-20 µL turquoise	manual single 2-20 µL yellow	manual single 10-100 µL	manual single 20-200 µL	manual single 100-1000 µL	manual multi 10 µL	manual multi 100 µL	manual multi 300 µL	electronic single 0.5-10 µL	electronic single 20-200 µL	electronic multi 300 µL
2-002-81-0	■	x						x			■		
2-003-C5-0	■	x						x			■		
2-062-81-0	■	x	■					■			■		
2-063-C5-0	■	x	■					■			■		
2-117-80-0		■	x	■									
2-118-C5-0		■	x	■									
2-119-80-0			x	■	■				■				
2-120-C5-0			x	■	■				■				
2-127-80-0			x	■	■	■						■	
2-128-C5-0			x	■	■	■						■	
2-134-80-0					■	■						■	■
2-135-C5-0					■	■						■	■
2-202-80-0							■						
2-203-C4-0							■						
2-204-80-0							■						
2-205-C4-0							■						
4-002-81-0	■	x						x			■		
4-003-C5-0	■	x						x			■		
4-062-81-0	■	x	■					■			■		
4-063-C5-0	■	x	■					■			■		
4-121-86-0		■	x	■	■	■			■			■	
4-122-C5-0		■	x	■	■	■			■			■	
4-127-80-0					■	■			■	■		■	■
4-128-C5-0					■	■			■	■		■	■
4-202-80-0							■						
4-203-C4-0							■						
4-204-80-0							■						
4-205-C4-0							■						
5-003-C5-0	■	x						x			■		
5-063-C5-0	■	x	■					■			■		
5-118-C5-0		■	x	■									
5-120-C5-0			x	■	■				■				
5-128-C5-0			x	■	■	■						■	
5-135-C5-0					■	■						■	■
5-203-C4-0							■						
5-205-C4-0							■						

Sartorius — Proline Plus / Picuss — manual single channel and electronic single and multi channel pipettes

Cat. No.	manual single 0.5-10 µL	manual single 2-20 µL	manual single 10-100 µL	manual single 20-200 µL	manual single 100-1000 µL	manual single 500-5000 µL	manual single 1000-10000 µL	electronic single 0.2-10 µL	electronic single 5-120 µL	electronic single 10-300 µL	electronic single 50-1000 µL	electronic multi 0.2-10 µL	electronic multi 5-120 µL	electronic multi 10-300 µL	electronic multi 50-1200 µL
2-002-81-0	■							■							
2-003-C5-0	■							■							
2-062-81-0	■							■				■			
2-063-C5-0	■							■				■			
2-117-80-0		■													
2-118-C5-0		■													
2-119-80-0			■												
2-120-C5-0			■												
2-127-80-0				■					x	■			x		
2-128-C5-0				■					x	■			x		
2-134-80-0				■										■	
2-135-C5-0				■										■	
2-202-80-0					■										
2-203-C4-0					■										
2-204-80-0					■										x
2-205-C4-0					■										x
4-002-81-0	■							■							
4-003-C5-0	■							■							
4-062-81-0	■							■				■			
4-063-C5-0	■							■				■			
4-121-86-0		■	■	■					■				x		
4-122-C5-0		■	■	■					■				x		
4-127-80-0		■	■	■						■				■	
4-128-C5-0		■	■	■						■				■	
4-202-80-0					■										
4-203-C4-0					■										
4-204-80-0					■										x
4-205-C4-0					■										x
5-003-C5-0	■							■							
5-063-C5-0	■							■				■			
5-118-C5-0		■													
5-120-C5-0			■												
5-128-C5-0				■					x	■			x		
5-135-C5-0				■										■	
5-203-C4-0					■										
5-205-C4-0					■										x

HTL — Labmate Pro (same as Labnet BioPette) — manual single and multi channel pipettes

Cat. No.	manual single 0.5-10 µL	manual single 2-20 µL	manual single 10-100 µL	manual single 20-200 µL	manual single 100-1000 µL	manual single 1000-5000 µL	manual single 1000-10000 µL	manual multi 0.5-10 µL	manual multi 50-300 µL
2-002-81-0	■							■	
2-003-C5-0	■							■	
2-062-81-0	■							■	
2-063-C5-0	■							■	
2-117-80-0		■							■
2-118-C5-0		■							■
2-119-80-0		■	■						■
2-120-C5-0		■	■						■
2-127-80-0		x	■	■					■
2-128-C5-0		x	■	■					■
2-134-80-0			■	■					■
2-135-C5-0			■	■					■
2-202-80-0					■				
2-203-C4-0					■				
2-204-80-0					■				
2-205-C4-0					■				
4-002-81-0	■							■	
4-003-C5-0	■							■	
4-062-81-0	■							■	
4-063-C5-0	■							■	
4-121-86-0		■	■	■					■
4-122-C5-0		■	■	■					■
4-127-80-0		■	■	■					■
4-128-C5-0		■	■	■					■
4-202-80-0					■				
4-203-C4-0					■				
4-204-80-0					■				
4-205-C4-0					■				
5-003-C5-0	■							■	
5-063-C5-0	■							■	
5-118-C5-0		■							■
5-120-C5-0		■	■						■
5-128-C5-0		x	■	■					■
5-135-C5-0			■	■					■
5-203-C4-0					■				
5-205-C4-0					■				

■ compatible x not compatible

+82-02-546-6211
+82-02-546-6212
www.bvskorea.com

ahn BVS